

Drug Test Series



Urine Drug Test

Premium Cup



■ MAX 22 test items

Select Cup



■ MAX 18 test items

Hex Cup




■ MAX 20 test items

Square Cup



■ MAX 18 test items



■ Multi-Panel Dipcard



■ Single Panel Dipcard

Oral Fluid Drug Test



D-Pen™ Mini



D-Pen™



Lollipop Test



Saltox

Hair Drug Test



Hair Drug Test

Powder/Liquid Drug Test



FTY Strip



XYL Strip



FTY + XYL Combo Strip

Available Test Items

Specimen	Test Item	Cut-off level	Format	Qualification
Urine	Amphetamine (AMP)	1000/500/300 ng/mL	Variable	Under 510k Application
	Barbiturates (BAR)	300 ng/mL	Variable	Under 510k Application
	Benzodiazepines (BZO)	300/200/100 ng/mL	Variable	Under 510k Application
	Buprenorphine (BUP)	10 ng/mL	Variable	Under 510k Application
	Cocaine (COC)	300/150/100 ng/mL	Variable	Under 510k Application
	Marijuana (THC)	50/40/25 ng/mL	Variable	Under 510k Application
	Methamphetamine (MET/mAMP)	1000/500/300 ng/mL	Variable	Under 510k Application
	Methadone (MTD)	300/200 ng/mL	Variable	Under 510k Application
	Methylenedioxymethamphetamine - ecstasy (MDMA)	500/300 ng/mL	Variable	Under 510k Application
	Morphine (MOP)	300/100 ng/mL	Variable	Under 510k Application
	Oxycodone (OXY)	100 ng/mL	Variable	Under 510k Application
	Phencyclidine (PCP)	25 ng/mL	Variable	Under 510k Application
	Propoxyphene (PPX)	300 ng/mL	Variable	Under 510k Application
	Tri-cyclic Antidepressants (TCA)	1000 ng/mL	Variable	Under 510k Application
	Synthetic Cannabis (K2)	50/25 ng/mL	Variable	FUO
	Opiate (OPI)	2000 ng/mL	Variable	Under 510k Application
	Tramadol HCl (TRA)	1000/200 ng/mL	Variable	FUO
	Methadone Metabolite (EDDP)	300 ng/mL	Variable	Under 510k Application
	Fentanyl (FTY)	10/20 ng/mL	Variable	FUO
	Cotinine (COT)	200 ng/mL	Variable	FUO
Ethyl Glucuronide (ETG)	500 ng/mL	Variable	FUO	
Ketamine (KET)	1000/500 ng/mL	Variable	FUO	
Monoacetylmorphine (6-MAM)	10 ng/mL	Variable	Under 510k Application	
Oral Fluid	Amphetamine (AMP)	50/40 ng/mL	Variable	EUO/FUO
	Barbiturates (BAR)	60/20 ng/mL	Variable	EUO/FUO
	Benzodiazepines (BZO)	30/10 ng/mL	Variable	FUO
	Buprenorphine (BUP)	5 ng/mL	Variable	FUO
	Cocaine (COC)	50/20 ng/mL	Variable	EUO/FUO
	Marijuana (THC)	50/40/12 ng/mL	Variable	EUO/FUO
	Methadone (MTD)	30 ng/mL	Variable	EUO/FUO
	Methamphetamine (MET/mAMP)	50 ng/mL	Variable	EUO/FUO
	Methylenedioxymethamphetamine (MDMA)	100/50/25 ng/mL	Variable	EUO/FUO
	Opiate (OPI)	50/40/15 ng/mL	Variable	EUO/FUO
	Morphine (MOP)	15 ng/mL	Variable	FUO
	Oxycodone (OXY)	20 ng/mL	Variable	EUO/FUO
	Phencyclidine (PCP)	10 ng/mL	Variable	EUO/FUO
	Cotinine (COT)	20/50 ng/mL	Variable	FUO
	Propoxyphene (PPX)	25 ng/mL	Variable	FUO
	Alcohol	0.02%	Variable	FUO
	Ethyl Glucuronide (ETG)	300/500 ng/mL	Variable	FUO
	Fentanyl (FTY)	20 ng/mL	Variable	FUO
	Monoacetylmorphine (6-MAM)	10 ng/mL	Variable	FUO
	Hair	Morphine (MOP)	0.2 ng/mg	Card
Methamphetamine (MET)		0.2 ng/mg	Card	FUO
Ketamine (KET)		0.2 ng/mg	Card	FUO
3 in 1 (MOP/MET/KET)		0.2 ng/mg	Card	FUO
Powder/Liquid	Fentanyl (FTY)	10 ng/mL	Strip	Harm Reduction
	Xylazine (XYL)	300 ng/mL	Strip	Harm Reduction
	FTY+ XYL combo test	10ng/mL + 300 ng/mL	Strip	Harm Reduction

Note: Additional tests arriving soon.